




CARE SURGICAL

CS Pad 2 inch




Place under patient supports to raise them up
when extra height is required.



Clean using standard
hospital wipes



Do not immerse in water



Do not autoclave



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, 12A Lower Main Street,
Lucan, Co. Dublin K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity
Assessed



CE Mark



5 060624 120893



LOT Number:

EC REP

Product Code: CSM-2568-2 Rev 2 01/24