



Technical Data Sheet

CSM-2646 CS Universal Arm Support



Description

- Reusable
- Attaches to most standard and specialist operating tables to support patients' arms in various positions during surgical procedures
- One-step set-up
- 3 rotating joints support full range of movement
- Central clamping mechanism for easy, fast and secure positioning of the arm support
- Including soft foam pad with padded strap
- Spherical joints for smooth rotation

Application Area

- OR
- Anaesthesia

Measurements

- L400 | W150 | H550 mm
- L15.7 | W5.9 | H21.7"

Material Composition

- Stainless steel
- PU-Foam (covered)
- Transfer coated PU with anti-bacterial treatment on knitted polyester

Cleaning

- Clean with standard hospital wipes

Order Information

Product Code	Product Name	Qty
CSM-2646	CS Universal Arm Support	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- Each unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- pressure ulcers
- neurological damage
- venous pressure, which can lead to ischaemia
- Do not use if there is any damage on the product
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624122125

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.02 3/24