



# Technical Data Sheet

## CSM2646 CS Universal Arm Support™



### Description

- Reusable
- Attaches to most standard and specialist operating tables to support patients' arms in various positions during surgical procedures
- One-step set-up
- 3 rotating joints support full range of movement
- Central clamping mechanism for easy, fast and secure positioning of the arm support
- Including soft foam pad with padded strap
- Spherical joints for smooth rotation
- Supports maximum weight of 10kg
- Product weight 4.3kg

### Application Area

- OR
- Anaesthesia

### Measurements

- Length: 460mm/18.11"
- Width: 150mm/5.9"
- Height: 550 mm/21.7"
- Max Height: (Rail to Top) 510mm/20.07"

### Material Composition

- Stainless steel
- PU-Foam (covered)
- Transfer coated PU with anti-bacterial treatment on knitted polyester

### Order Information

Product Code	Product Name	Qty
CSM2646	CS Universal Arm Support™	1 Unit

### Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- Each unit is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

### Contraindications

Incorrect positioning of the patient on this product can cause:

- pressure ulcers
- neurological damage
- venous pressure, which can lead to ischaemia
- Do not use if there is any damage on the product
- Do not autoclave
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624122125

### Cleaning

- Clean with standard hospital wipes

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson  
MA 01749, USA  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark