




CARE SURGICAL

CS Supine Top -3" Foam pad



CS Supine Top; reusable covered foam support.

Reusable foam pad for pressure management, for use with CS Supine Top- Base Only or as replacement pads for compatible Allen or Jackson style supine base*.

*Allen® is the registered trademark of Allen Medical Systems Inc.



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
Marlborough
MA 01752
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, Office 2, 12a Lower Main Street,
Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



UK Conformity
Assessed



Reuse



CE Mark



Latex Free



Non Sterile



5 060624 122194

LOT Number:




Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia

Product Code: CSM-2722-3 Rev.2 11 /23




CARE SURGICAL


CS Supine Top -3" Foam pad



- 1** Place the attachment on the frame so that the brackets are over the bars.



- 2** Snap the locking clip shut to secure the attachment in position.



- 3** Place the CS Supine Top Foam Pad in the centre of the CS Supine Top Base.


Cleaning instructions



Clean using standard hospital wipes



Do not immerse in water



Do not autoclave



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
Marlborough
MA 01752
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, Office 2, 12a Lower Main Street,
Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



UK Conformity
Assessed



Reuse



CE Mark



Latex Free



Non Sterile



5 060624 122194

EC REP

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia

Product Code: CSM-2722-3 Rev.2 11 /23