




CARE SURGICAL
Specialists in Patient Positioning

**NEW
PRODUCT**

CS Supine Top (Base Only)

The CS Supine Top (Base Only) is a solid base that can be attached to most standard dual-rail operating tables* on the market.


- + Can be used e.g. with the different CS Supine Top foam pads available in 1, 2 or 3" thickness
- + With the help of the Supine Top Base, dual rail tables can be easily converted into tables with a solid surface
- + The cut-outs enable additional attachment options for e.g. table rail adapters or other accessories for more flexibility



CS Supine Top (Base Only)


Code: CSM-2722-7
Qty: 1

Related products:



**CS Supine Top- 1''
Foam Pad**

Code: CSM-2722-1
Qty: 1 Unit




**CS Supine Top- 2''
Foam Pad**

Code: CSM-2722-2
Qty: 1 Unit



**CS Supine Top- 3''
Foam Pad**

Code: CSM-2722-3
Qty: 1 Unit



**CS Supine Top-
with 2'' Foam Pad**

Code: CSM-2722
Qty: 1 Set

*Compatible with rail sizes that measure 1.25" (3.2cm) wide and 1.5" (3.8 cm) tall with 14.5" (36.8 cm) space between.



Contact us for more information:
www.care-surgical.com | info@care-surgical.com