



Technical Data Sheet

CSM2722-7 CS Supine Top - Base Only



Description

- Sectional mattress base for use on dual rail operating table
- Reusable
- Compatible with the Allen style dual-rail operating tables
- Sectional base to support pressure relieving pads for patients in the supine or prone position.
- Clean using non-alcohol hospital wipes
- Number of uses: Unrestricted, subject to damage (see damage indicators in the Contraindications column)

Application Area

- OR
- Anaesthesia

Measurements

- Length: 500mm/19.69"
- Width: 305mm/12.01"
- Height: 75mm/2.95"

Material Composition

- Nylon
- Stainless steel
- ABS Plastic
- Carbone Fibre
- Latex & phthalate free

Cleaning

- Clean with standard hospital wipes

Order Information

Product Code	Product Name	Qty
CSM2722-7	CS Supine Top - Base Only	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight

Contraindications

- Ensure the platform is securely attached to the table before use
- Check the platform for damage before and after each use
- Do not use if damaged, contact supplier for advice
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624122392

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Australian Sponsor:

MDSS AU Pty Ltd, Ground Floor, Suite 55,
97-99 Bathurst Street, Sydney, NSW 2000, Australia

Ref: CSM2722-7- TDS (Rev.04) 03/2026