



# Technical Data Sheet

## CSM2722-2 CS Supine Top - 2" Foam Pad



### Order Information

| Product Code | Product Name                | Qty    |
|--------------|-----------------------------|--------|
| CSM2722-2    | CS Supine Top - 2" Foam Pad | 1 Unit |

### Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 1 unit is packaged into a cardboard box for sale

### Description

- Sectional pressure management pads for patients in the prone, supine or lateral position
- Compatible with the Allen-Style Dual-Rail operating tables
- For use on the CSM2722-7 CS Supine Top- Base Only
- Number of uses: Unrestricted, subject to damage (see damage indicators in the Contraindications column)

### Application Area

- OR
- Anaesthesia

### Measurements

- Length: 460mm/18.11"
- Width: 305mm/12.01"
- Height: 60mm/2.36"

### Material Composition

- Polyurethane foam
- Latex & phthalate free
- Transfer coated PU with anti-bacterial treatment on knitted polyester

### Cleaning

- Clean using non-alcohol hospital wipes

### Storage Instructions

- First in first out
- Store in a cool dry place
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

### Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia

Examine the pad for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.

- Do not autoclave
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624122187

### Labelling




Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson  
MA 01749, USA  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor:

MDSS AU Pty Ltd, Ground Floor, Suite 55,  
97-99 Bathurst Street, Sydney, NSW 2000, Australia



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark