



Technical Data Sheet

CSM-2666 CS Traction Boot Disposable



Description

- Traction boot insert made from polyurethane foam
- Single use - disposal after use with regular clinical waste
- Cushioned insert compatible with the Allen Traction Boot
- Aids patient comfort
- Reduces the incidence of pressure related injury whilst the patient's foot and lower leg are supported in the Allen Traction Boot
- Allen® is the registered trademark of Allen Medical Systems Inc

Application Area

- OR
- Anaesthesia

Measurements

- Length 660mm/25.98"
- Width 610mm/24.02"
- Height 8mm/0.31"

Material Composition

- Polyurethane foam
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM-2666	CS Traction Boot Disposable	6 Units

Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 6 units are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia
- Reuse of this product could lead to cross infection
- Do not use if product is torn or damaged
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120985

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023