



# Technical Data Sheet

## CSM-2666 CS Traction Boot Disposable



### Description

- Traction boot insert made from polyurethane foam
- Single use - disposal after use with regular clinical waste
- Cushioned insert compatible with the Allen Traction Boot
- Aids patient comfort
- Reduces the incidence of pressure related injury whilst the patient's foot and lower leg are supported in the Allen Traction Boot
- Allen® is the registered trademark of Allen Medical Systems Inc

### Application Area

- OR
- Anaesthesia

### Measurements

- Length 660mm/25.98"
- Width 610mm/24.02"
- Height 8mm/0.31"

### Material Composition

- Polyurethane foam
- Hook & Loop
- Latex & phthalate free

### Order Information

| Product Code | Product Name                | Qty     |
|--------------|-----------------------------|---------|
| CSM-2666     | CS Traction Boot Disposable | 6 Units |

### Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 6 units are packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

### Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia
- Reuse of this product could lead to cross infection
- Do not use if product is torn or damaged
- Do not autoclave
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120985

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
Marlborough  
MA 01752  
United States  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical  
PO Box 60, Landsborough  
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023