



Technical Data Sheet

CSM-2667 CS Head Strap Disposable



Description

- Secures patient's head and chin on beach chair type tables
- Head strap made from polyurethane foam laminated with brushed nylon and with an added simple hook and loop fixing
- Single use - disposal after use with regular clinical waste
- Straps are designed to distribute force and help reduce the prevalence of pressure injuries

Application Area

- OR
- Anaesthesia

Measurements

- Length 870mm/34.25"
- Width 120mm/4.72"
- Height 20mm/0.79"

Material Composition

- Brushed nylon
- Velcro fixing
- Polyurethane foam
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM-2667	CS Head Strap Disposable	10 Sets

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 10 sets packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- Incorrect positioning of the patient on this product can cause pressure marking or ulcers.
- Securing this product too loosely could result in movement of the head during the procedure
- Do not use if this product is torn or damaged
- Reuse of this product could lead to cross infection
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120992

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main
Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not
Reuse



Latex Free



Non Sterile



UK Conformity
Assessed



CE Mark

Rev.01 11/09/2023