



Technical Data Sheet

CSM-2722 CS Supine Top - With 2" Foam Pad



Description

- Reusable
- Sectional pressure management pads for patients in the supine or prone position.
- Compatible with the Allen-Style Dual-Rail operating tables.
- For use on the Supine base plate.
- Clean using non-alcohol hospital wipes
- Number of uses: Unrestricted, subject to damage (see damage indicators in the Contraindications column)
- Compatible with rail sizes that measure 1.25" (3.2cm) wide and 1.5" (3.8 cm) tall

Application Area

- OR
- Anaesthesia

Measurements

Base (CSM-2722-7) -

- Length: 500mm/19.69"
- Width: 305mm/12.01"
- Height: 75mm/ 2.95"

Top (CSM-2722-2) -

- Length: 460mm/18.11"
- Width: 305mm/12.01"
- Height: 60mm/ 2.36"

Material Composition

- Polyurethane foam
- Latex & phthalate free

Cleaning

- Clean with standard hospital wipes

Order Information

Product Code	Product Name	Qty
CSM-2722	CS Supine Top - With 2" Foam Pad	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia

Examine the pad for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.

- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624127229

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023