



Technical Data Sheet

CSM-2588 CS Memory Foam Sheet (200cm x 55cm x 5cm)



Order Information

Product Code	Product Name	Qty
CSM-2588	CS Memory Foam Sheet (200cm x 55cm x 5cm)	2 Units

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 2 units are packaged into a cardboard box for sale

Description

- For use on operating tables
- Can be cut to size
- Reduces and distributes pressure to reduce the risk of pressure ulcers
- The undulating structure of the foam enables air to circulate to assist in maintaining the patient's temperature

Application Area

- OR
- Anaesthesia

Measurements

- Length 2000mm/79"
- Width 550mm/22"
- Height 50mm/2"

Material Composition

- Memory foam
- Latex & phthalate free

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624122552

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023