



Technical Data Sheet

CSM2739-1 CS Table Rail Adaptor EU



Description

- Adaptor for use on standard dual-rail tables, with EU table rail. It is used to provide a secure method of mounting (holding) other medical devices, such as armboards
- Reusable - clean using alcohol-based hospital wipes
- For use with dual rail tables, compatible with CS Carbon Spine Extension Frame

Application Area

- OR
- Anaesthesia

Measurements

- Length 155mm/6.10"
- Width 85mm/3.35"
- Height 70mm/2.76"

Material Composition

- Aluminium,
- Stainless steel
- Acetal
- Nylon
- Latex & phthalate free

Cleaning

- Clean using alcohol-based hospital wipes

Order Information

Product Code	Product Name	Qty
CSM2739-1	CS Table Rail Adaptor EU	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place

Contraindications

- Incorrect attachment to dual rail tables or the CS Carbon Spine Extension Frame could result in the attachments failing
- Do not use if the adaptor appears to be defective or damaged
- Only to be used by trained medical personnel
- Do not exceed the device safe working load
- Do not immerse in water
- Do not autoclave
- After each use, clean the device with alcohol-based wipes
- Make sure that the device is dry before you store it or use it again

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN: 5060624123573

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark