



CS Standard Replacement Stirrup Pads

Replacement stirrup pads for Allen Medical OEM Stirrups. Secure and cushion the patient's feet and legs during lithotomy positioning, reducing the chance of pressure related injury. Maximum patient weight of 350lbs (159kg).

- + Made from supportive memory foam, fully encapsulating the patient's foot
- + Compatible with Allen® Medical OEM Stirrups
- + Reusable, easy wipe-clean

CS Standard Replacement Pads

Code: CSM-2750

Qty: 1 Pair



Contact us for more information:
www.care-surgical.com | info@care-surgical.com