



Technical Data Sheet

CSM2750 CS Standard Replacement Stirrup Pads (Pair)



Description

- Secures and cushions patient's feet and legs in lithotomy stirrups
- Compatible with the Allen® Medical OEM Stirrups
- For use in a variety of surgical procedures
- Reduces the incidence of pressure ulcers and marking
- Made from supportive memory foam, fully encapsulating the patient's foot
- Maximum patient weight: 350lbs (159kg)

Application Area

- OR
- Anaesthesia

Measurements

- Length 570mm/22.4"
- Width 560mm/22"
- Height 100mm/3.9"

Material Composition

- Polyurethane foam
- Latex & phthalate free
- Transfer coated PU with anti-bacterial treatment on knitted polyester

Cleaning

- Clean with non-alcohol hospital wipes

Order Information

Product Code	Product Name	Qty
CSM2750	CS Standard Replacement Stirrup Pads (Pair)	1 Pair

Packaging Information

- Each Pair is packaged in a recyclable polyethene bag
- 1 pair is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can lead to temporary or permanent paralysis
- Venous pressure, which can lead to ischaemia.
- Examine the product for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1 Basic
- UDI/GTIN No: 5060624123689

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity
Assessed



CE Mark