



# Technical Data Sheet

## CSM-2763 CS Din Rail Clamp - EU



### Description

- Clamp with standard rail for attachment to standard dual rail table – perfect for attaching arm boards, retractors and other similar accessories
- Easy to attach and fix with a screw
- The body of the clamp fits perfectly into the cut out of the CS Supine Top Base CSM-2722-7
- Compatible with rail sizes that measure 1.25" (3.2cm) wide and 1.5" (3.8 cm) tall
- Weight capacity: 250kg

### Application Area

- OR
- Anaesthesia

### Measurements

- Length 135mm/5"
- Width 160mm/6"
- Height 115mm/4.5"

### Material Composition

- Aluminium
- Stainless Steel
- Plastic

### Order Information

Product Code	Product Name	Qty
CSM-2763	CS Din Rail Clamp - EU	1 Unit

### Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit are packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight

### Contraindications

- Incorrect attachment to the CS Carbon Spine Extension could result in the attachments failing
- Do not use if the adaptor appears to be defective or damaged.
- Only to be used by trained medical personnel
- Do not exceed the device safe working load
- Do not immerse in water
- Do not autoclave
- After each use, clean the device with alcohol-based wipes.
- Make sure that the device is dry before you store it or use it again.

### Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624124068
- Latex & phthalate free

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
Marlborough  
MA 01752  
United States  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 10/07/2024