



# CS Din Rail Clamp

Clamp with din rail (EU) for attachment to standard dual rail table – perfect for attaching arm boards, retractors and other similar accessories.

- + Easy to attach and fix with a screw
- + The body of the clamp fits perfectly into the cut out of the CS Supine Top Base CSM-2722-7
- + Compatible with rail sizes that measure 1.25" (3.2cm) wide and 1.5" (3.8 cm) tall

## CS Din Rail Clamp - EU

Code: CSM-2763

Qty: 1 Unit



## Related Products



CS Supine Top -  
Base Only

Code: CSM-2722-7

Qty: 1 Unit



CS Universal Arm  
support - Moulded

Code: CSM-2647

Qty: 1 Unit



**Contact us for more information:**  
[www.care-surgical.com](http://www.care-surgical.com) | [info@care-surgical.com](mailto:info@care-surgical.com)