




# CARE SURGICAL

## CS Table Rail Adaptor-EU



- 1** Position the rail adaptor in the desired position.
- 2** Place the rail adaptor flush to the operating table rail. Ensure the adaptor is placed with the clip on the inside facing the operating table and the adaptor facing outwards.
- 3** Close the secure locking clip to secure the adaptor ready for use.



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
Tel: 01704 336671

USA: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson,  
MA 01749, USA  
Tel: 866-243-4107



Maximum  
Weight



Medical Device



UK Conformity  
Assessed



Reuse



CE Mark



Latex Free



Non Sterile



5 060624 123573



LOT Number:

**www.care-surgical.com**

Care Surgical Supplies Limited, Office 2, 12a Lower Main  
Street, Lucan Co. Dublin, K78 X5P8 Ireland



Product Code: CSM-2739-1 Rev 4 01/24