



CARE SURGICAL
Specialists in Patient Positioning

CS Lateral Draw Sheet

A contoured sheet made from polyurethane foam, layered with soft brushed nylon, and equipped with a reliable Velcro fastening. Ensures smooth transitions from supine to lateral to prone positions on the Allen lateral top.

- + Controlled turn of the patient from supine to lateral to prone position
- + Latex & phthalate free
- + Single use to aid infection control

More Products from Our Disposable Range:



CS Adjustable Prone Leg Cushion

Code: CSM-2590
Qty: 2 Sets



CS Traction Boot Foam Insert

Code: CSM-2666
Qty: 6 Units




CS Disposable Leg Holder

Code: CSM-2738
Qty: 12 Units



CS Prone Chest and Pelvic Rolls

Code: CSM-2902
Qty: 2 Units




CS Standard Arm Cushions

Code: CSM-2575
Qty: 10 Units



CS ProneOne Head Support

Code: CSM-2520
Qty: 8 Units



CS Lateral Draw Sheet

Code: CSM-2672
Qty: 12 Units

Allen® is the registered trademark of Allen Medical Systems Inc.

Contact us for more information:
www.care-surgical.com | info@care-surgical.com