




# CARE SURGICAL

# CS Disposable Leg Holder


CS Disposable Leg Holder - For use with compatible arthroscopic leg holder to provide safe and comfortable support.



**1** Position the leg holder pad under the patients leg with the white side facing outwards.



**2** Wrap the leg holder pad securely around the patients leg and secure using the strap.



**3** Use with compatible arthroscopic leg holder.



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
Tel: 01704 336671

USA: Care Surgical LLC  
Marlborough  
MA 01752  
Tel: 866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)

Care Surgical Supplies Limited, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



UK Conformity Assessed



Do Not Reuse



CE Mark



Latex Free



Non Sterile



5 060624 127380



Australian Sponsor, Summit Surgical  
PO Box 60, Landsborough  
QLD 4550 Australia

LOT Number:

Product Code: CSM-2738 Rev 1 09/22