




Technical Data Sheet

CSM2674 CS Prone Arm Secure - Large



Order Information

Product Code	Product Name	Qty
CSM2674	CS Prone Arm Secure - Large	10 Units

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 10 units are packaged into a cardboard box for sale

Description

- To secure the patient's arms by their side in the prone position
- Sheet made from brushed nylon on 2mm polyurethane foam with velcro fixings
- Wraps under and over the patient's arms and secured with Velcro to hold their arms securely by their side in the prone position

Application Area

- OR
- Anaesthesia

Measurements

- Length: 1350mm/53"
- Width: 450mm/18"
- Height: 5mm/0.2"

Material Composition

- Polyurethane foam
- Brushed nylon
- Velcro fixings
- Latex & phthalate free

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight

Contraindications

- Incorrect positioning of the arms can cause damage to the radial or ulner nerve
- Securing the sheet too tightly and cause discomfort or tissue damage to the patient
- Securing the sheet too loosely could cause the patient's to move during the procedure
- Do not use if this product is torn or damaged
- Reuse of this product could lead to cross infection
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624122071

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Ref: CSM2674 - TDS (Rev.03) 03/2026