




Technical Data Sheet

CSM2738 CS Disposable Leg Holder



Description

- Leg holder made from polyurethane foam with added hook and loop fixing
- For use with compatible Arthroscopic Leg Holder Systems during knee arthroscopy procedures
- Cushions, secures and stabilises the leg
- Single use - disposal after use with regular clinical waste

Application Area

- OR
- Anaesthesia

Measurements

- Length: 960mm/37.80"
- Width: 90mm/3.54"
- Height: 10mm/0.39"

Material Composition

- Polyurethane foam
- Velcro fixing
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM2738	CS Disposable Leg Holder	10 Sets

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 10 sets are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- Incorrect positioning of the patient in this product can cause pressure marking or ulcers
- Securing this product too loosely could result in movement of the leg during the procedure
- Do not use if this product is torn or damaged
- Reuse of this product could lead to cross infection
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN: 5060624127380

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark