




# Technical Data Sheet

## CSM-2674 CS Prone Arm Secure - Large



### Description

- To secure the patient's arms by their side in the prone position
- Sheet made from brushed nylon on 2mm polyurethane foam with velcro fixings
- Wraps under and over the patient's arms and secured with Velcro to hold their arms securely by their side in the prone position

### Application Area

- OR
- Anaesthesia

### Measurements

- Length 1350mm/53"
- Width 450mm/18"
- Height 5mm/0.2"

### Material Composition

- Polyurethane foam
- Brushed nylon
- Velcro fixings
- Latex & phthalate free

### Order Information

Product Code	Product Name	Qty
CSM-2674	CS Prone Arm Secure - Large	10 Units

### Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 10 units are packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight

### Contraindications

- Incorrect positioning of the arms can cause damage to the radial or ulnar nerve
- Securing the sheet too tightly and cause discomfort or tissue damage to the patient
- Securing the sheet too loosely could cause the patient's to move during the procedure
- Do not use if this product is torn or damaged
- Reuse of this product could lead to cross infection
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624122071

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
Marlborough  
MA 01752  
United States  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main  
Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical  
PO Box 60, Landsborough  
QLD 4550 Australia



Medical Device



Do Not  
Reuse



Latex Free



Non Sterile



UK Conformity  
Assessed



CE Mark

Rev.01 11/09/2023