




Technical Data Sheet

CSM-2738 CS Disposable Leg Holder



Description

- Leg holder made from polyurethane foam with added hook and loop fixing
- For use with compatible Arthroscopic Leg Holder Systems during knee arthroscopy procedures
- Cushions, secures and stabilises the leg
- Single use - disposal after use with regular clinical waste

Application Area

- OR
- Anaesthesia

Measurements

- Length 960mm/37.80"
- Width 90mm/3.54"
- Height 10mm/0.39"

Material Composition

- Polyurethane foam
- Velcro fixing
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM-2738	CS Disposable Leg Holder	10 Sets

Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 10 sets are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- Incorrect positioning of the patient in this product can cause pressure marking or ulcers
- Securing this product too loosely could result in movement of the leg during the procedure
- Do not use if this product is torn or damaged
- Reuse of this product could lead to cross infection
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624127380

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023