



Technical Data Sheet

TDS CSM-2568-1, -2, -3 CS Pads



Description

- Reusable
- Placed under patient supports to raise them up when extra height is required
- To be used e.g. in combination with CS Prone Head Support or CS ProneOne

Application Area

- OR
- Anaesthesia

Measurements

- -1: L275 | W275 | H25 mm / L10.8 | W10.8 | H1"
- -2: L275 | W275 | H50 mm / L10.8 | W10.8 | H2"
- -3: L275 | W275 | H76 mm / L10.8 | W10.8 | H3"

Material Composition

- PU-Foam (covered)

Cleaning

- Clean with standard hospital wipes

Order Information

Product Code	Product Name	Qty
CSM-2568-1	CS Pad 1"	1 Unit
CSM-2568-2	CS Pad 2"	1 Unit
CSM-2568-3	CS Pad 3"	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- Each unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- pressure ulcers
 - neurological damage
 - venous pressure, which can lead to ischaemia
-
- Do not use if there is any damage on the product
 - Do not autoclave
 - Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No - 1": 5060624122262
- Basic UDI/GTIN No - 2": 5060624120893
- Basic UDI/GTIN No - 3": 5060624122279

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity
Assessed



CE Mark

Rev.01 8/23