



# Technical Data Sheet

## CSM-2568-3 CS Pad 3"



### Description

- Flat pad made from polyurethane foam with a wipeable skin covering
- Raises the height of patient head supports to aid in achieving a neutral neckline
- Placed under patient supports to raise them up when extra height is required
- Number of uses: Unrestricted, subject to damage (see damage indicators in the Contraindications column)

### Application Area

- OR
- Anaesthesia

### Measurements

- Length 300mm/11.81"
- Width 275mm/10.83"
- Height 76mm/2.99"

### Material Composition

- Polyurethane foam
- Latex & phthalate free
- Wipeable skin covering
- Transfer coated PU with anti-bacterial treatment on knitted polyester

### Cleaning

- Clean with standard hospital wipes

### Order Information

Product Code	Product Name	Qty
CSM-2568-3	CS Pad 3"	1 Unit

### Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 1 unit is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

### Contraindications

- The pad and support may become unstable if placed on an uneven surface
- Not recommended for use with supports larger than 300x300mm
- Examine the pad for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624122279

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
Marlborough  
MA 01752  
United States  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main  
Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity  
Assessed



CE Mark

Rev.02 04/03/2024