




CARE SURGICAL


Product Code: CSM2568-3

CS Pad 3 inch

⚠️ Only to be used by trained clinicians. Ensure product is secure before positioning patient.




Cleaning instructions:



Clean using standard hospital wipes



Do not immerse in water



Do not autoclave

Place under patient supports to raise them up when extra height is required.

⚠️ Do NOT put patient directly onto pad



UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107



www.care-surgical.com

Care Surgical Supplies Limited, 12A Lower Main Street,
Lucan, Co. Dublin K78 X5P8 Ireland



⚠️ **SAFE DISPOSAL:** customers should adhere to any local laws and regulations to the safe disposal of medical devices.
⚠️ **DO NOT USE** products that appear to be damaged or worn. In an event of product damage or patient injury, please contact Care Surgical.



Ref: CSM2568-3 - IFU (Rev.04) 03/2025