




CARE SURGICAL
Specialists in Patient Positioning

The CS Carbon Spine Extension - Schmitz™

The CS Carbon Spine Extension - Schmitz™ can be fitted to the Schmitz Diamond Operating Table offering a patient position and imaging capabilities comparable to that of a specialist dual-rail spine table.



Find out more:





CARE SURGICAL
Specialists in Patient Positioning

Surgeries requiring a break - Microdiscectomy and Osteotomy



Table can be flexed and reflexed allowing the intervertebral space to be opened and closed with the hand control.



The CS Carbon Spine Extension is fully compatible with Wilson style frames



CARE SURGICAL
Specialists in Patient Positioning

Anterior and Lateral Approach



The CS Carbon Fibre Supine Tops allow the patient to be positioned supine for anterior approach surgery such as ALIF



The CS Carbon Fibre Supine Tops allow the patient to be positioned laterally enabling lateral approach surgery such as XLIF



The CS Carbon Spine Extension can also be utilised for pain management along with other procedures which require an imaging table allowing for 360 degree artifact free imaging.



Whats Included:



CS Carbon Spine Extension - Schmitz™ (Set)

Code: CSM2932
Qty: 1 Set

Includes:

CSM-2730 CS Carbon Spine Extension Frame
CSM-2737 CS Carbon Spine Extension Leg
CSM-2953 CS Carbon Spine Extension Adaptor - Schmitz



CS Prone Plus Head Support System™

Code: CSM2565
Qty: 1 System
(Helmet & Mirror)



CS Universal Arm Support - Moulded™

Code: CSM2647
Qty: 1 Unit



CS Easy Shin Support™

Code: CSM3020
Qty: 1 Unit



CS Knee Rings (Pair)

Code: CSM2821
Qty: 1 Pair



CS Table Rail Adaptor - EU

Code: CSM2739-1
Qty: 1 Unit



CS Prone Plus Face Cushion

Code: CSM2512
Qty: 10 Units



CS Comfort Covers - 5 Piece

Code: CSM2600
Qty: 6 Units



CS Head Tray Attachment - Lowered

Code: CSM2751
Qty: 1 Unit



CS Chest Support

Code: CSM2709
Qty: 1 Unit



CS Hip Supports (pair)

Code: CSM2717
Qty: 1 Pair



CS Popper Hips (Pair)

Code: CSM2743
Qty: 1 Set



CS Thigh Supports (pair)

Code: CSM2718
Qty: 1 Pair



CS Popper Thigh Supports (pair)

Code: CSM2744
Qty: 1 Pair



CS Pad 2"

Code: CSM2568-2
Qty: 1 Unit

Optional Accessories:



CS Din Rail Clamp

Code: CSM2763
Qty: 1 Unit



CS Supine Top-Base Only

Code: CSM2722-7
Qty: 1 Unit



CS Supine Top - 3" Foam Pad

Code: CSM2568-2
Qty: 1 Unit



CS Comfort Covers - 5 Piece

Code: CSM2600
Qty: 6 Units



CS Carbon Spine Extension Cart

Code: 3084
Qty: 1 Unit

Contact us for more information:
www.care-surgical.com | info@care-surgical.com