



Technical Data Sheet

CSM2724 CS Stirrup Cart



Description

- The Stirrup Leg holders can be easily hooked into the handle of the cart
- The clamps can be stowed in the lower part of the cart
- Four antistatic double castors allow easy transport
- All castors can be braked for optimal safety
- Easy to clean with wipe disinfectant

Application Area

- OR
- Anaesthesia

Measurements

- Length 570mm/22.44"
- Width 480mm/18.90"
- Height 1050mm/41.34"

Material Composition

- Latex & phthalate free
- Plastic
- Stainless Steel
- Polypropylene

Order Information

Product Code	Product Name	Qty
CSM2724	CS Stirrup Cart	1 Unit

Packaging Information

- Each cart is wrapped with protective foam and packaged in a recyclable polyethene bag
- Each cart is packaged into a cardboard box for sale.

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Transport and store in dry conditions

Contraindications

- Not loaded with more than 30kg
- Always push, not pull to avoid tipping over
- Always hook the stirrups properly to prevent them from falling out
- Always brake the castors when the cart is not in motion to avoid rolling away and accidents
- Do not damage surfaces with sharp objects

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624123672

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.02 01/04/2025