



CS Blade Clamp & Blade Clamp Plus

Reusable clamp for securing blade-style accessories to standard operating tables via the table rails.

- + Used to secure surgical table accessories with flat posts to table rails
- + Easily secured with screw system for stable fixing of operating accessories
- + Reusable - clean using alcohol-based hospital wipes
- + **CS Blade Clamp:** Slides on for easy fitting to standard operating tables
- + **CS Blade Clamp Plus:** Can be attached and locked at any position on the table rail



CS Blade Clamp Plus- US

Code: CSM-2762
Qty: 1 Unit



CS Blade Clamp Plus- EU

Code: CSM-2767
Qty: 1 Unit



CS Blade Clamp - UK/EU

Code: CSM-2761
Qty: 1 Unit



Contact us for more information:
www.care-surgical.com | info@care-surgical.com