



CARE SURGICAL

Specialists in Patient Positioning

CS Easy Prone

The CS Easy Prone simplifies the positioning of the chest and pelvis in the prone position.

- + Fully radiolucent and fully sealed
- + The special shape prevents pressure in the armpit, neck and genital area
- + Enables compression-free positioning of the abdomen
- + Suitable for many different body shapes, as it can be flexibly adjusted
- + Allows movements of the table without compromises
- + Enables a natural lordosis of the spine
- + Available in various heights for optimal positioning



CS Easy Prone - 150mm - Chest & Pelvis Set

Code: CSM-2699
Qty: 1 Set

CS Easy Prone - 170mm - Chest & Pelvis Set

Code: CSM-2918
Qty: 1 Set


CS Easy Prone - 200mm - Chest & Pelvis Set

Code: CSM-2919
Qty: 1 Set



CS Prone Plus Head Support System

Code: CSM-2565
Qty: 1 System
(Helmet and Mirror)



Use with the CS Prone Plus Head Support System

Contact us for more information:
www.care-surgical.com | info@care-surgical.com