



**CARE SURGICAL**  
Specialists in Patient Positioning

## CS Head Tray Attachment - Lowered

The CS Head Tray Attachment - Lowered is compatible with most standard dual-rail operating systems on the market.

- + Maximum flexibility in height adjustment thanks to the lowered design
- + Quick-release clipping mechanism for easy attachment and adjustment



\*Compatible with rail sizes that measure 1.25" (3.2cm) wide and 1.5" (3.8 cm) tall with 14.5" (36.8cm) space between.

### CS Head Tray Attachment

Code: CSM-2751  
Qty: 1 Unit

## Related Products:



### CS Prone Plus Head Support System

Code: CSM-2565  
Qty: 1 System  
(Helmet & Mirror)



### CS Prone Head Support System

Code: CSM-2560  
Qty: 1 System  
(Helmet & Mirror)




### CS Prone Plus Face Cushion

Code: CSM-2512  
Qty: 10 Units



### CS Prone Small Face Cushion

Code: CSM-2514  
Qty: 10 Units



### CS Pad 2"

Also available in 1" & 3"

Code: CSM-2568-2  
Qty: 1 Unit



**Contact us for more information:**  
[www.care-surgical.com](http://www.care-surgical.com) | [info@care-surgical.com](mailto:info@care-surgical.com)