



CARE SURGICAL

Specialists in Patient Positioning

INSTRUCTIONS FOR USE

CSM2643

CS Supine & Lateral Head Support




CARE SURGICAL

Product Code: CSM2643

CS Supine & Lateral Head Support

⚠ Only to be used by trained clinicians. Ensure product is secure before positioning patient.



The CS Supine & Lateral Head Support provides safe and comfortable support to patients heads in the supine and lateral position.

Cleaning instructions:



Clean using standard hospital wipes



Do not immerse in water



Do not autoclave

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

Australian Sponsor:

MDSS AU Pty Ltd
Ground Floor, Suite 55, 97-99 Bathurst Street,
Sydney, NSW 2000, Australia

www.care-surgical.com

Care Surgical Supplies Limited, Office 2, 12a Lower Main Street,
Lucan Co. Dublin, K78 X5P8 Ireland



SAFE DISPOSAL: customers should adhere to any local laws and regulations to the safe disposal of medical devices.



DO NOT USE products that appear to be damaged or worn. In an event of product damage or patient injury, please contact Care Surgical.



EU REP




CARE SURGICAL

Artikelnummer: CSM2643

CS Supine & Lateral Head Support

⚠ Nur für den Einsatz durch geschultes Personal geeignet. Vergewissern Sie sich, dass das Produkt sicher befestigt ist, bevor Sie den Patienten positionieren.



Das CS Supine & Lateral Head Support bietet den Patienten in Rücken- und Seitenlage eine sichere und bequeme Positionierung des Kopfes.

Reinigungshinweise:



Mit handelsüblichen
Krankenhaus-Reinigungstüchern
reinigen.



Nicht in Wasser eintauchen.



Nicht autoklavieren.



UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107



www.care-surgical.com

Care Surgical Supplies Limited, Office 2, 12a Lower Main Street,
Lucan Co. Dublin, K78 X5P8 Ireland



Sichere Entsorgung: Kunden sollten alle lokalen Gesetze und Vorschriften zur sicheren Entsorgung von Medizinprodukten einhalten.



Verwenden Sie keine Produkte: die beschädigt oder abgenutzt erscheinen. Im Falle einer Produktbeschädigung oder einer Verletzung des Patienten wenden Sie sich bitte an Care Surgical.





CARE SURGICAL
Specialists in Patient Positioning