



Technical Data Sheet

CSM2751 CS Head Tray Attachment - Lowered



Order Information

Product Code	Product Name	Qty
CSM2751	CS Head Tray Attachment - Lowered	1 Unit

Description

- To be used with a head support and/or face cushion to support a patient's head
- Secure the tray to the side bars of the operating table and lock onto place with the quick-release locking clips
- Reusable - clean using non-alcohol hospital wipes

Application Area

- OR
- Anaesthesia

Measurements

- Length: 501mm/20"
- Width: 296mm/12"
- Height: 155mm/6"

Material Composition

- Aluminium
- Stainless Steel
- Carbon Fibre plates
- Latex and phthalate free

Cleaning

- Clean using non-alcohol hospital wipes

Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- Ensure the platform is securely attached to the table before use
- Check the platform for damage before and after each use
- Do not use if damaged Contact supplier for advice
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN: 5060624127519

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main
Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity
Assessed



CE Mark