



Technical Data Sheet

CSD-150-S CS Torso Support Cushion - Small



Description

- To support the patient's torso in the prone position
- Single use - disposal after use with regular clinical waste
- Raises the patient's torso in the prone position
- To assist in venous return and diaphragmatic excursion return
- Achieving a natural neutral neckline
- For patients weighing up to 200kg

Application Area

- OR
- Anaesthesia

Measurements

- Length 440mm/17.32"
- Width 340mm/13.38"
- Height 140mm/5.51"

Material Composition

- Polyurethane foam
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSD-150-S	CS Torso Support Cushion - Small	2 Units

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 2 units are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

- Incorrect positioning of the patient can cause:
- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia;
- Venous pressure, which can lead to ischaemia
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624122675

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023