



Technical Data Sheet

CSM-2589 CS Tourniquet Pad Disposable - 900mm x 130mm x 6mm



Description

- Padding to help protect patient's arms from pinching and shear forces
- Single Use
- Protects patient's arms from pinching and shear forces from the use of a tourniquet

Application Area

- OR
- Anaesthesia

Measurements

- Length 900mm/35.43"
- Width 130mm/5.12"
- Height 6mm/0.24"

Material Composition

- Polyurethane foam
- Latex & phthalate free

Order Information

Product Code	Product Name	Qty
CSM-2589	CS Tourniquet Pad Disposable - 900mm x 130mm x 6mm	100 Units

Packaging Information

- Each bag of 10 units are packaged in a recyclable polyethylene bag
- 10 bags are packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624122569

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 08/09/2023