



Technical Data Sheet

CSM-2643 CS Supine & Lateral Head Support



*Also available in Small. Image for size comparison.

Description

- Supports patients' heads in the supine or lateral position
- Reduces the incidence of pressure ulcers and marking
- Reusable

Application Area

- OR
- Anaesthesia

Measurements

- Length 231mm/9.09"
- Width 267mm/10.51"
- Height 95mm/3.74"

Material Composition

- Polyurethane foam
- Latex & phthalate free

Cleaning

- Clean using alcohol-based hospital wipes

Order Information

Product Code	Product Name	Qty
CSM-2643	CS Supine & Lateral Head Support	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Do not store anything heavy on top
- Avoid contact with sharp objects

Contraindications

- Incorrect positioning of the patient on the CS Supine & Lateral Head Support can cause:
- Pressure marking or ulcers
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel
- Do not immerse in water

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624120688

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity
Assessed



CE Mark

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main
Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia

Rev.01 11/09/2023