




# Technical Data Sheet

## CSM-2644 CS Supine & Lateral Head Support - Small



### Order Information

Product Code	Product Name	Qty
CSM-2644	CS Supine & Lateral Head Support - Small	1 Unit

\*Also available in Large. Image for size comparison.

### Description

- Supports patients' heads in the supine or lateral position
- Reduces the incidence of pressure ulcers and marking
- Reusable
- Suitable for smaller and / or paediatric patients

### Application Area

- OR
- Anaesthesia

### Measurements

- Length 157.5mm/6.20"
- Width 180.5mm/7.11"
- Height 62mm/2.44"

### Material Composition

- Polyurethane foam
- Latex & phthalate free

### Cleaning

- Clean using alcohol-based hospital wipes

### Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Do not store anything heavy on top
- Avoid contact with sharp objects

### Contraindications

- Incorrect positioning of the patient on the CS Supine & Lateral Head Support - Small can cause:
- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia;
- Venous pressure, which can lead to ischaemia.
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel
- Do not immerse in water

### Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624122644

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
Marlborough  
MA 01752  
United States  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical  
PO Box 60, Landsborough  
QLD 4550 Australia



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 23/11/2023