



# Technical Data Sheet

## CSM-2733-5 CS Armboard - 3"/80mm Pad - with strap



### Description

- For use on standard compatible armboards to support patients' arms
- Reusable
- Suitable to support patients' arms in a wide range of positions
- Padded strap with hook and loop fixings wraps around the patient's arm to hold it in position during the procedure and then around the armboard to secure the pad in place.
- Reusable: Clean using non-alcohol hospital wipes
- Number of uses: Unrestricted, subject to damage (see damage indicators in the Contraindications column)

### Application Area

- OR
- Anaesthesia

### Measurements

- Length 670mm/17.32"
- Width 140mm/5.51"
- Height 80mm/0.79"
- Strap: 500mm/ 19.68"

### Material Composition

- Latex & phthalate free
- Transfer coated PU with anti-bacterial treatment on knitted polyester

### Cleaning

- Clean with standard hospital wipes

### Order Information

Product Code	Product Name	Qty
CSM-2733-5	CS Armboard - 3"/80mm Pad - with strap	1 Unit

### Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

### Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia

Examine the pad for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.

- Do not autoclave
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624122118

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
Marlborough  
MA 01752  
United States  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.02 04/03/2024