



# Technical Data Sheet

## CSM-2741 CS Ultrafin compatible Replacement Pads (pair)



### Description

- Reusable
- To protect patient's legs in stirrups
- Secures and cushions patient's feet and legs in lithotomy stirrups
- Compatible with the Ultrafin stirrups
- For use in a variety of surgical procedures
- Reduces the incidence of pressure ulcers and marking
- Easy wipe-clean cover
- Maximum patient weight: 350lbs (159kg)

### Application Area

- OR
- Anaesthesia

### Measurements

- Length 725mm/29"
- Width 620mm/24"
- Height 30mm/1.2"

### Material Composition

- Reusable wipeable cover
- Hook & Loop fixings
- Polyurethane foam core
- Latex & phthalate free
- Transfer coated PU with anti-bacterial treatment on knitted polyester

### Cleaning

- Clean using non-alcohol hospital wipes

### Order Information

Product Code	Product Name	Qty
CSM-2741	CS Ultrafin compatible Replacement Pads (pair)	1 Pair

### Packaging Information

- Each pair is packaged in a recyclable polyethene bag
- 1 pair is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

### Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia
- Venous pressure, which can lead to ischaemia

Examine the pad for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced

- Do not autoclave
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624122149

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
Marlborough  
MA 01752  
United States  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical  
PO Box 60, Landsborough  
QLD 4550 Australia



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 11/09/2023