




Technical Data Sheet

CSM-2907 CS Universal Half Roll



Description

- Prone or Supine Positioning support
- Supine: Distributes pressure above the hollow of the patients knee and helps to prevent the knee from being pushed through. Reduces the occurrence of pressure related injuries
- Prone: Helps to position the patients feet without any pressure on the arch of the foot and toe
- Comfortable for both conscious and unconscious patients
- Reusable

Application Area

- OR
- Anaesthesia

Measurements

- Length 500mm/19.69"
- Width 180mm/7.09"
- Height 120mm/4.72"

Material Composition

- Polyurethane foam
- Latex & phthalate free
- PU-Foam cover - wipeable

Cleaning

- Clean using alcohol-based hospital wipes

Order Information

Product Code	Product Name	Qty
CSM-2907	CS Universal Half Roll	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

- Incorrect positioning of the patient on this product could result in the patient experiencing:
- Pressure ulcers, neurological damage, which can lead to temporary or permanent paralysis and venous pressure, which can lead to ischaemia.
- Examine the product for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624129070

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical
PO Box 60, Landsborough
QLD 4550 Australia



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 14/12/2023