



# Technical Data Sheet

## CSM-2919 CS Easy Prone - 200mm - Chest & Pelvis Support



### Order Information

Product Code	Product Name	Qty
CSM-2919	CS Easy Prone - 200mm - Chest & Pelvis Support	1 Set

### Description

- Set consists of: CSM-2924 CS Easy Prone - 200mm - Chest Only & CSM-2925 CS Easy Prone - 200mm - Pelvis Only
- Raises the patient's chest and pelvis in the prone position
- Enables compression-free positioning of the abdomen
- Reusable, fully radiolucent and fully sealed
- The design prevents pressure in the armpit, neck and genital area.
- Suitable for many different body shapes, as it can be flexibly adjusted
- Enables a natural lordosis of the spine
- Allows movements of the table without compromises

### Application Area

- OR
- Anaesthesia

### Measurements

#### Chest:

- Length 290mm/11.42"
- Width 500mm/19.69"
- Height 200mm/7.87"

#### Pelvis:

- Length 376mm/14.8"
- Width 500mm/19.69"
- Height 200mm/7.87"

### Material Composition

- Polyurethane foam
- Latex & phthalate free
- Weft knitted polyamide, coated with anti-microbial polyurethane coating

### Packaging Information

- Each component of the set is packaged in a recyclable polyethylene bag
- 1 set is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place

### Contraindications

- Incorrect positioning of the patient on the CS Easy Prone can cause:
- Pressure marking or ulcers
- Neurological damage, which can cause paraplegia;
- Venous pressure, which can lead to ischaemia.
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624129193

### Cleaning

- Clean using alcohol-based hospital wipes

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
Marlborough  
MA 01752  
United States  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)

EC REP

Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.01 23/11/2023