



# Technical Data Sheet

## CSM2644 CS Supine & Lateral Head Support - Small



\*Also available in Large. Image for size comparison.

### Description

- Supports patients' heads in the supine or lateral position
- Reduces the incidence of pressure ulcers and marking
- Reusable
- Suitable for smaller and / or paediatric patients

### Application Area

- OR
- Anaesthesia

### Measurements

- Length: 157.5mm/6.20"
- Width: 180.5mm/7.11"
- Height: 62mm/2.44"

### Material Composition

- Polyurethane foam
- Latex & phthalate free

### Cleaning

- Clean using alcohol-based hospital wipes

### Order Information

Product Code	Product Name	Qty
CSM2644	CS Supine & Lateral Head Support - Small	1 Unit

### Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Do not store anything heavy on top
- Avoid contact with sharp objects

### Contraindications

- Incorrect positioning of the patient on the CS Supine & Lateral Head Support - Small can cause:
  - Pressure marking or ulcers
  - Neurological damage, which can cause paraplegia;
  - Venous pressure, which can lead to ischaemia.
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel
- Do not immerse in water

### Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624122644

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson  
MA 01749, USA  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark