



Technical Data Sheet

CSM2907 CS Universal Half Roll



Description

- Ideal for positioning patients in the supine position
- Suitable for a wide range of patients
- Reduces the occurrence of pressure-related injuries
- Reusable

Application Area

- OR
- Anaesthesia

Measurements

- Length: 500mm/19.69"
- Width: 150mm/7.09"
- Height: 75mm/4.72"

Material Composition

- Polyurethane foam
- Latex & phthalate free
- PU-Foam cover - wipeable

Cleaning

- Clean using non-alcohol hospital wipes

Order Information

Product Code	Product Name	Qty
CSM2907	CS Universal Half Roll	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethylene bag
- 1 unit is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

- Incorrect positioning of the patient on this product could result in the patient experiencing:
- Pressure ulcers, neurological damage, which can lead to temporary or permanent paralysis and venous pressure, which can lead to ischaemia.
- Examine the product for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN: 5060624129070

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
Marlborough
MA 01752
United States
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main
Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity
Assessed



CE Mark

Ref:CSM2907 – TDS (Rev.03) 03/2026