




CARE SURGICAL

Product Code: CSM2644

CS Supine & Lateral Head Support - Small

 Only to be used by trained clinicians. Ensure product is secure before positioning patient.



The CS Supine & Lateral Head Support – Small provides safe and comfortable support for the heads of smaller patients (e.g. paediatric patients) in both the supine and lateral positions.

Cleaning instructions:



Clean using standard hospital wipes



Do not immerse in water



Do not autoclave



UK: Care Surgical Ltd
Unit 6 Ringtail Road,
Burscough,
L40 8JY
Tel: 01704 336671



USA: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107



www.care-surgical.com

Care Surgical Supplies Limited, Office 2, 12a Lower Main Street,
Lucan Co. Dublin, K78 X5P8 Ireland



 **SAFE DISPOSAL:** customers should adhere to any local laws and regulations to the safe disposal of medical devices.
 **DO NOT USE** products that appear to be damaged or worn. In an event of product damage or patient injury, please contact Care Surgical.



5 060624 122644

Ref: CSM2644 - IFU (Rev.05) 03/2025