



# Technical Data Sheet

## CSM2589 CS Tourniquet Pad Disposable - 900mm x 130mm x 6mm



### Description

- Padding to help protect patient's arms from pinching and shear forces
- Single Use
- Protects patient's arms from pinching and shear forces from the use of a tourniquet

### Application Area

- OR
- Anaesthesia

### Measurements

- Length 900mm/35.43"
- Width 130mm/5.12"
- Height 6mm/0.24"

### Material Composition

- Polyurethane foam
- Latex & phthalate free

### Order Information

Product Code	Product Name	Qty
CSM2589	CS Tourniquet Pad Disposable - 900mm x 130mm x 6mm	100 Units

### Packaging Information

- Each bag of 10 units are packaged in a recyclable polyethene bag
- 10 bags are packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

### Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Do not use if this product is torn or damaged
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- Basic UDI/GTIN No: 5060624122569

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson  
MA 01749, USA  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland

Australian Sponsor, Summit Surgical  
PO Box 60, Landsborough  
QLD 4550 Australia



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Rev.02 01/04/2025