



CARE SURGICAL

Specialists in Patient Positioning

INSTRUCTIONS FOR USE

CSM3020 **CS Easy Shin Support™**




CARE SURGICAL

Product Code: CSM3020

CS Easy Shin Support™

⚠ Only to be used by trained clinicians. Ensure product is secure before positioning patient.



CS Easy Shin Support™

Exposed dorsum of the foot and protection of the dorsalis pedis artery.

Ensure feet do not touch operating table.

Cleaning instructions:



Clean using standard hospital wipes



Do not immerse in water



Do not autoclave



UK: Care Surgical Ltd
Unit 6 Ringtail Road,
Burscough,
L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, Office 2, 12a Lower Main Street,
Lucan Co. Dublin, K78 X5P8 Ireland



SAFE DISPOSAL: customers should adhere to any local laws and regulations to the safe disposal of medical devices.
DO NOT USE products that appear to be damaged or worn.
In an event of product damage or patient injury, please contact Care Surgical.



EU REP




CARE SURGICAL

Artikelnummer: CSM3020

CS Easy Shin Support™

⚠ Nur für den Einsatz durch geschultes Personal geeignet. Vergewissern Sie sich, dass das Produkt richtig liegt, bevor Sie den Patienten positionieren.



CS Easy Shin Support™

Gewährleistet freiliegende Fußrück-
en sowie Schutz der Arteria dorsalis
pedis. Achten Sie darauf, dass die
Füße den Operationstisch nicht
berühren.


Reinigungshinweise:



Mit handelsüblichen
Krankenhaus-Reinigungs-
tüchern reinigen



Nicht in Wasser
eintauchen



Nicht autoklavieren



UK: Care Surgical Ltd
Unit 6 Ringtail Road,
Burscough,
L40 8JY
Tel: 01704 336671

USA: Care Surgical LLC
43 Broad Street
Suite B106, Hudson,
MA 01749, USA
Tel: 866-243-4107

www.care-surgical.com

Care Surgical Supplies Limited, Office 2, 12a Lower Main Street,
Lucan Co. Dublin, K78 X5P8 Ireland



SICHERE ENTSORGUNG: Kunden sollten alle lokalen Gesetze und
Vorschriften zur sicheren Entsorgung von Medizinprodukten einhalten.
VERWENDEN SIE KEINE PRODUKTE die beschädigt oder abgenutzt
erscheinen. Im Falle einer Produktbeschädigung oder einer Verletzung
des Patienten wenden Sie sich bitte an Care Surgical.



EU REP



CARE SURGICAL
Specialists in Patient Positioning