



Technical Data Sheet

CSM3024 CS Universal Supine Head Support



Order Information

Product Code	Product Name	Qty
CSM3024	CS Universal Supine Head Support	1 Unit

Packaging Information

- Each unit is packaged in a recyclable polyethene bag
- 1 unit are packaged into a cardboard box for sale

Description

- Reusable - clean using non-alcohol hospital wipes
- Secures the patients head in supine, lithotomy or trendelenburg position
- Reduces the incidence of pressure ulcers and marking
- For use in a variety of surgical procedures
- Compatible for example with the CSM2648 CS Supine Bariatric Support

Application Area

- OR
- Anaesthesia

Measurements

- Length: 240mm/9.45"
- Width: 210mm/8.27"
- Height: 110mm/4.33"

Material Composition

- Polyurethane foam
- Latex & phthalate free
- Transfer coated PU with anti-bacterial treatment on knitted polyester

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Do not use if this product is torn or damaged
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/TGA/UKCA:Y
- MD Class: 1
- UDI-DI/GTIN: 5060624125812

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark