



Technical Data Sheet

CSM3024 CS Universal Supine Head Support



Description

- Reusable - clean using non-alcohol hospital wipes
- Secures the patients head in supine, lithotomy or trendelenburg position
- Reduces the incidence of pressure ulcers and marking
- For use in a variety of surgical procedures
- Compatible for example with the CSM2648 CS Supine Bariatric Support

Application Area

- OR
- Anaesthesia

Measurements

- Length: 240mm/9.45"
- Width: 210mm/8.27"
- Height: 110mm/4.33"

Material Composition

- Polyurethane foam
- Latex & phthalate free
- Transfer coated PU with anti-bacterial treatment on knitted polyester

Cleaning

- Clean with standard hospital wipes

Order Information

Product Code	Product Name	Qty
CSM3024	CS Universal Supine Head Support	1 Unit

Packaging Information

- Each Pair is packaged in a recyclable polyethene bag
- 1 pair is packaged into a cardboard box for sale

Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top
- Avoid storing where sharp objects could pierce the skin of the product

Contraindications

Incorrect positioning of the patient on this product can cause:

- Pressure marking or ulcers
- Neurological damage, which can lead to temporary or permanent paralysis
- Venous pressure, which can lead to ischaemia. Examine the product for damage before and after each use. If the skin is pierced or torn, it cannot be cleaned effectively and should be replaced.
- Do not immerse in water
- Do not autoclave
- Only to be used by trained medical personnel

Further Information

- Disposable/Reusable: Reusable
- Sterility: N
- CE/FDA/TGA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN: 5060624125812

Labelling



Manufacturer

UK: Care Surgical Ltd
Unit 6 Ringtail Road
Burscough, L40 8JY
United Kingdom
Tel: +44 1704 336671

US: Care Surgical LLC
43 Broad Street
Suite B106, Hudson
MA 01749, USA
Tel: +1-866-243-4107

www.care-surgical.com



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark