



# Technical Data Sheet

## CSM2926 CS Easy Prone Comfort Covers



### Description

- Disposable covers for the CS Easy Prone (all sizes)
- Helps protect the CS Easy Prone from damage or contamination
- Latex and phthalate free
- Single-use aids infection control

### Application Area

- OR
- Anaesthesia

### Measurements

#### Chest

- Length: 500mm/19.69"
- Width: 250mm/9.84"
- Height: 170mm/6.69"

#### Pelvis

- Length: 500mm/19.69"
- Width: 335mm/13.19"
- Height: 170mm/6.69"

### Material Composition

- White foam: Polyurethane
- Blue non-woven (waterproof): 55% Woodpulp  
45% Polyester
- Latex & phthalate free

### Order Information

Product Code	Product Name	Qty
CSM2926	CS Easy Prone Comfort Covers	6 Sets

### Packaging Information

- Each set is packaged in a recyclable polyethene bag
- 6 sets is packaged into a cardboard box for sale

### Storage Instructions

- First in first out
- Store in a cool dry place
- Avoid direct sunlight
- Do not store anything heavy on top

### Contraindications

- If incorrectly applied to the CS Easy Prone, this product can cause pressure marks or ulcers.
- Do not reuse
- Only to be used by trained medical personnel

### Further Information

- Disposable/Reusable: Disposable
- Sterility: N
- CE/FDA/UKCA: Y
- MD Class: 1
- UDI-DI/GTIN No: 5060624129261

### Labelling



Manufacturer

UK: Care Surgical Ltd  
Unit 6 Ringtail Road  
Burscough, L40 8JY  
United Kingdom  
Tel: +44 1704 336671

US: Care Surgical LLC  
43 Broad Street  
Suite B106, Hudson  
MA 01749, USA  
Tel: +1-866-243-4107

[www.care-surgical.com](http://www.care-surgical.com)



Care Surgical Supplies Ltd, Office 2, 12a Lower Main Street, Lucan Co. Dublin, K78 X5P8 Ireland



Medical Device



Do Not Reuse



Latex Free



Non Sterile



UK Conformity Assessed



CE Mark

Ref: CSM2926 - TDS (Rev.04) 03/2026